


GSK3117391

Product Code: F009092
 CAS No.: 1018673-42-1
 Mol. Weight: C22H33N3O4
 Mol. Formula: 403.52
 Storage Temp: Powder: -20°C, 36 months; 4°C, 24 months
 In solvent: -80°C, 24 months; -20°C, 12 months
 Synonyms: ESM-HDAC391; CHR-5154; HDAC-IN-3
 Target: HDAC
 Pathway: Cell Cycle/DNA Damage; Epigenetics



Biological Activity

制备储备液:

溶剂体积 浓度 \ 质量	1 mg	5 mg	10 mg
1 mM	2.4782 mL	12.3910 mL	24.7819 mL
5 mM	0.4956 mL	2.4782 mL	4.9564 mL
10 mM	0.2478 mL	1.2391 mL	2.4782 mL

Bioactive: HDAC-IN-3 is a histone deacetylase (HDAC) inhibitor, extracted from patent WO/2008040934 A1. Target: HDAC



访问 FELIX 中国官网获取最新优惠活动、产品动态
 FELIX 所有产品只能用于科学研究
 任何个人和单位不能将产品用于人体、临床等非科学研究领域和用途