


(±)-2-Methylbutyryl chloride-d9

Product Code: F078513
 CAS No.: 1219795-10-4
 Mol. Weight: C5D9C10
 Mol. Formula: 129.63
 Storage Temp: Please store the product under the recommended conditions in the Certificate of Analysis.
 Target: Isotope-Labeled Compounds
 Pathway: Others



Biological Activity

制备储备液:

溶剂体积 浓度	质量 1 mg	5 mg	10 mg
1 mM	7.7143 mL	38.5713 mL	77.1426 mL
5 mM	1.5429 mL	7.7143 mL	15.4285 mL
10 mM	0.7714 mL	3.8571 mL	7.7143 mL

Bioactive: (±)-2-Methylbutyryl chloride-d9 is the deuterium labeled (±)-2-Methylbutyryl chloride.

In Vitro: Stable heavy isotopes of hydrogen, carbon, and other elements have been incorporated into drug molecules, largely as tracers for quantitation during the drug development process. Deuteration has gained attention because of its potential to affect the pharmacokinetic and metabolic profiles of drugs.



访问 FELIX 中国官网获取最新优惠活动、产品动态
FELIX 所有产品只能用于科学研究
任何个人和单位不能将产品用于人体、临床等非科学研究领域和用途